

FTT

Four-Way Float Top Table

This sturdy table is ideally suited for a variety of general radiographic needs including small hospitals, clinics, ambulatory care centers and private offices. The floating table top glides freely on a precision roller bearing system and can support patients up to 500 lbs. (227 kg).



*Foot treadle electric lock
release runs the entire length
of pedestal.*

Product Features

- 17" x 17" image receptor with heavy-duty manual tray
- Heavy-duty electric locks
- All steel structure and welded base provides maximum strength and reliability
- Smooth, stain resistant fiber resin top
- Pedestal base can house generator electronics to save space

